

For Research Use Only

VSIG1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87189-2-PBS



Basic Information

Catalog Number:

87189-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6740

GenBank Accession Number:

BC043216

GeneID (NCBI):

340547

UNIPROT ID:

Q86XK7

Full Name:

V-set and immunoglobulin domain containing 1

Calculated MW:

387 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

252228B2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02900-1, Human VSIG1 Recombinant Matched Antibody Pair - PBS only. 87189-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6740. 87189-1-PBS was HRP conjugated as the detection antibody. Range: 6.25-400 ng/mL.