

For Research Use Only

SLC12A1/NKCC2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87210-1-PBS



Basic Information

Catalog Number:

87210-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5325

GenBank Accession Number:

BC040138

GeneID (NCBI):

6557

UNIPROT ID:

Q13621

Full Name:

solute carrier family 12
(sodium/potassium/chloride
transporters), member 1

Calculated MW:

121 kDa

Observed MW:

140-180 kDa

Purification Method:

Protein A purification

CloneNo.:

252363D10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated and 37°C incubated mouse kidney tissue were subjected to SDS PAGE followed by western blot with 87210-1-RR (SLC12A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87210-1-PBS in a different storage buffer formulation.