

For Research Use Only

MYO18A Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87218-3-PBS



Basic Information

Catalog Number:

87218-3-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC039612

GeneID (NCBI):

399687

UNIPROT ID:

Q92614

Full Name:

myosin XVIII A

Calculated MW:

233 kDa

Purification Method:

Protein A purification

CloneNo.:

251878D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MYO18A concentration was determined to be 3.4 ng/mL in HaCaT cell extract based on a 1.4 mg/mL extract load, HL-60 cell extract were served as a negative control.



Cytometric bead array standard curve of MP02925-1, MYO18A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87218-3-PBS. Detection antibody: 87218-2-PBS. Standard: DNA011-LNP. Range: 0.156-20 ng/mL