

For Research Use Only

MYO18A Recombinant monoclonal antibody, PBS Only

Catalog Number:87218-4-PBS



Basic Information

Catalog Number:

87218-4-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC039612

GeneID (NCBI):

399687

UNIPROT ID:

Q92614

Full Name:

myosin XVIII A

Calculated MW:

233 kDa

Observed MW:

230-245 kDa

Purification Method:

Protein A purification

CloneNo.:

251878A9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87218-4-RR (MYO18A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87218-4-PBS in a different storage buffer formulation.