

For Research Use Only

BRCA2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87225-2-PBS



Basic Information

Catalog Number:

87225-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31064

GenBank Accession Number:

NM_000059

GeneID (NCBI):

675

UNIPROT ID:

P51587

Full Name:

breast cancer 2, early onset

Calculated MW:

384 kDa

Purification Method:

Protein A purification

CloneNo.:

252307C7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean BRCA2 concentration was determined to be 4.91 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load and 2.22 ng/mL in A2780 cell extract based on a 1.20 mg/mL extract load.



Plasma of twenty-three individual healthy human donors was measured. The BRCA2 concentration of detected samples was determined to be 0.77 ng/mL with a range of 0.05-6.07 ng/mL.



Sandwich ELISA standard curve of MP03287-1, Human BRCA2 Recombinant Matched Antibody Pair - PBS only. 87225-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag31064. 87225-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.