

For Research Use Only

PDIA6 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87234-2-PBS



Basic Information

Catalog Number:

87234-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5454

GenBank Accession Number:

NM_005742.3

GeneID (NCBI):

10130

UNIPROT ID:

Q15084-1

Full Name:

protein disulfide isomerase family A,
member 6

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

252495B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PDIA6 concentration was determined to be 2378.7 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1285.7 ng/mL in HCT 116 cell extract based on a 1.3 mg/mL extract load, 257.2 ng/mL in hTERT-RPE1 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02930-1, PDIA6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87234-3-PBS. Detection antibody: 87234-2-PBS. Standard: Eg5454. Range: 0.625-80 ng/mL.