

For Research Use Only

# PDIA6 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87234-3-PBS



## Basic Information

Catalog Number:

87234-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5454

GenBank Accession Number:

NM\_005742.3

GeneID (NCBI):

10130

UNIPROT ID:

Q15084-1

Full Name:

protein disulfide isomerase family A, member 6

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

252495E12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PDIA6 concentration was determined to be 2378.7 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1285.7 ng/mL in HCT 116 cell extract based on a 1.3 mg/mL extract load, 257.2 ng/mL in hTERT-RPE1 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02930-1, PDIA6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87234-3-PBS. Detection antibody: 87234-2-PBS. Standard: Eg5454. Range: 0.625-80 ng/mL.