

For Research Use Only

MUL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87235-1-PBS



Basic Information

Catalog Number:

87235-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9119

GenBank Accession Number:

BC010101

GeneID (NCBI):

79594

UNIPROT ID:

Q969V5

Full Name:

mitochondrial E3 ubiquitin ligase 1

Calculated MW:

352 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

252482F6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MUL1 concentration was determined to be 893.20 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load, 617.00 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 477.35 ng/mL in THP-1 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02942-1, Human MUL1 Recombinant Matched Antibody Pair - PBS only. 87235-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag9119. 87235-1-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL