

For Research Use Only

XRCC1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87251-1-PBS



Basic Information

Catalog Number: 87251-1-PBS	GenBank Accession Number: BC023593	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7515	CloneNo.: 252518B12
Isotype: IgG	UNIPROT ID: P18887	
Immunogen Catalog Number: AG15172	Full Name: X-ray repair complementing defective repair in Chinese hamster cells 1	
	Calculated MW: 633 aa, 70 kDa	
	Observed MW: 80-90 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

XRCC1 (X-ray Repair Cross Complementing group 1) is a molecular scaffold protein that interacts with multiple enzymatic components of DNA single-strand break repair (SSBR) including DNA kinase, DNA phosphatase, DNA polymerase, DNA deadenylase, and DNA ligase activities that collectively are capable of accelerating the repair of a broad range of DNA single-strand breaks (SSBs) (PMID: 18515191). XRCC1 is 633 amino acids in length which is comprised of an N-terminal domain (NTD) of ~160 amino acids that interacts with Pol β , a central BRCT domain of ~90 amino acids that interacts with PARP1 (PMID: 9584196). This antibody is not suitable for murine samples.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87251-1-RR (XRC1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87251-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using XRC1 antibody (87251-1-RR, Clone: 252518B12) at dilution of 1:800 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 87251-1-PBS in a different storage buffer formulation.



Bi-layer interferometry (BLI) kinetic assays of 87251-1-RR against Human XRC1 were performed. The affinity constant is below 1 pM.