

For Research Use Only

GP1BA Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87285-5-PBS



Basic Information

Catalog Number:

87285-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3051

GenBank Accession Number:

NM_000173.7

GeneID (NCBI):

2811

ENSEMBL Gene ID:

ENSG00000185245

UNIPROT ID:

P07359

Full Name:

glycoprotein Ib (platelet), alpha polypeptide

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

252200B8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of two individual healthy human donors was measured. The CD42b concentrations in the detected samples were 290.39 ng/mL and 184.18 ng/mL, respectively.



Sandwich ELISA standard curve of MP02948-2, Human CD42b Recombinant Matched Antibody Pair - PBS only. 87285-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3051. 87285-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.