

For Research Use Only

PALMD Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87300-2-PBS



Basic Information

Catalog Number: 87300-2-PBS	GenBank Accession Number: BC012570	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 54873	CloneNo.: 252551F6
Isotype: IgG	UNIPROT ID: Q9NP74	
Immunogen Catalog Number: AG9796	Full Name: palmdelphin	
	Calculated MW: 551 aa, 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PALMD concentration was determined to be 148.3 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 42.8 ng/mL in Mouse lung tissue extract based on a 5.0 mg/mL extract load.



Cytometric bead array standard curve of MP02954-1, PALMD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87300-2-PBS. Detection antibody: 87300-1-PBS. Standard: Ag9796. Range: 0.469-60 ng/mL.