

For Research Use Only

Mouse CD316/IGSF8 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87310-2-PBS



Basic Information

Catalog Number: 87310-2-PBS	GenBank Accession Number: AF439263.1	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 140559	CloneNo.: 252500E12
Isotype: IgG	UNIPROT ID: Q8R366	
Immunogen Catalog Number: EG6320	Full Name: immunoglobulin superfamily, member 8	
	Calculated MW: 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IGSF8/CD316 concentration was determined to be 865.1 ng/mL in mouse brain tissue extract based on a 2.2 mg/mL extract load, A431 cell extract were served as a negative control.



Cytometric bead array standard curve of MP02958-1, MOUSE IGSF8/CD316 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87310-1-PBS. Detection antibody: 87310-2-PBS. Standard: Eg6320. Range: 0.625-80 ng/mL.