

For Research Use Only

INPP5F Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87311-2-PBS



Basic Information

Catalog Number: 87311-2-PBS	GenBank Accession Number: BC111493	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 22876	CloneNo.: 252661A12
Isotype: IgG	UNIPROT ID: Q9Y2H2	
Immunogen Catalog Number: AG39933	Full Name: inositol polyphosphate-5-phosphatase F	
	Calculated MW: 1132 aa, 128 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean INPP5F concentration was determined to be 59.4 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02978-1, INPP5F Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87311-3-PBS. Detection antibody: 87311-2-PBS. Standard: Ag39933. Range: 0.938-60 ng/mL.