

For Research Use Only

IFNW1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87318-2-PBS



Basic Information

Catalog Number: 87318-2-PBS	GenBank Accession Number: BC117292	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3467	CloneNo.: 251487A9
Isotype: IgG	UNIPROT ID: P05000	
Immunogen Catalog Number: EG5386	Full Name: interferon, omega 1	
	Calculated MW: 22 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of sixteen individual healthy human donors was measured. The IFNWI1 concentration of detected samples was determined to be 24.1 pg/mL with a range of 4.4-56.7 pg/mL.



Sandwich ELISA standard curve of MP03285-1, Human IFNWI1 Recombinant Matched Antibody Pair - PBS only. 87318-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5386. 87318-2-PBS was HRP conjugated as the detection antibody. Range: 31.3-1000 pg/mL.