

For Research Use Only

PTPRJ/CD148 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87321-3-PBS



Basic Information

Catalog Number:

87321-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3492

GenBank Accession Number:

NM_002843.4

GeneID (NCBI):

5795

UNIPROT ID:

Q12913-1

Full Name:

protein tyrosine phosphatase, receptor type, J

Calculated MW:

146 kDa

Purification Method:

Protein A purification

CloneNo.:

252382E12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PTPRJ/CD148 concentration was determined to be 61.1 ng/mL in HeLa cell extract based on a 2.2 mg/mL extract load, A-204 cell extract were served as a negative control.



Cytometric bead array standard curve of MP02965-1, PTPRJ/CD148 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87321-3-PBS. Detection antibody: 87321-2-PBS. Standard: Eg3492. Range: 0.781-100 ng/mL