

For Research Use Only

TLE1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87332-5-PBS



Basic Information

Catalog Number:

87332-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17299

GenBank Accession Number:

BC010100

GeneID (NCBI):

7088

UNIPROT ID:

Q04724

Full Name:

transducin-like enhancer of split 1
(E(sp1) homolog, Drosophila)

Calculated MW:

770 aa, 83 kDa

Purification Method:

Protein A purification

CloneNo.:

252574A6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TLE1 concentration was determined to be 4.12 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 1.89 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load and 1.59 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02977-2, Human TLE1 Recombinant Matched Antibody Pair - PBS only. 87332-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17299. 87332-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL