

For Research Use Only

MYOTI Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87334-2-PBS



Basic Information

Catalog Number:

87334-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6068

GenBank Accession Number:

/

GeneID (NCBI):

9499

UNIPROT ID:

Q9UBF9

Full Name:

myotilin

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

252646B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MYOT concentration was determined to be 89.1 ng/mL in Mouse skeletal muscle tissue extract based on a 4.1 mg/mL extract load.



Cytometric bead array standard curve of MP02976-1, MYOT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87334-1-PBS. Detection antibody: 87334-2-PBS. Standard: Eg6068. Range: 0.781-100 ng/mL.