

For Research Use Only

Melan-A Recombinant monoclonal antibody, PBS Only

Catalog Number: 87351-2-PBS



Basic Information

Catalog Number: 87351-2-PBS	GenBank Accession Number: BC014423	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2315	CloneNo.: 252647G1
Isotype: IgG	UNIPROT ID: Q16655	
Immunogen Catalog Number: AG13346	Full Name: melan-A	
	Calculated MW: 13 kDa	
	Observed MW: 13-20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Melan-A is a melanocyte differentiation antigen, recognized by autologous cytotoxic T lymphocytes. Melan-A is also called MART-1 (melanoma antigen recognized by T cells). The Melan-A/MART-1 gene encodes this protein, 20-22 kDa, associated with endoplasmic reticulum and melanosomes. The function of the protein is unknown. Melan-A, isolated as a melanoma-specific antigen, is a transmembrane protein, which is expressed in skin, retina and the majority of cultured melanocytes as well as in melanomas. Melan A is expressed in more than 85% of melanomas.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SK-MEL-30 cells were subjected to SDS PAGE followed by western blot with 87351-2-RR (MLANA antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87351-2-PBS in a different storage buffer formulation.