

For Research Use Only

# Mouse HGFR/c-Met Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87352-2-PBS



## Basic Information

<b>Catalog Number:</b> 87352-2-PBS	<b>GenBank Accession Number:</b> /	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 17295	<b>CloneNo.:</b> 252395A7
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P16056	
<b>Immunogen Catalog Number:</b> EG6351	<b>Full Name:</b> met proto-oncogene	
	<b>Calculated MW:</b> 154 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




Serum of six mice was measured. The mouse HGFR/c-Met concentration of detected samples was determined to be 90.31 ng/mL with a range of 75.06 - 109.96 ng/mL.



Sandwich ELISA standard curve of MP02974-1, Mouse HGFR/c-Met Recombinant Matched Antibody Pair - PBS only. 87352-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6351. 87352-1-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL



The mean HGFR/c-Met concentration was determined to be 86.1 ng/mL in Mouse serum.



Cytometric bead array standard curve of MP02974-1, MOUSE HGFR/c-Met Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87352-2-PBS. Detection antibody: 87352-1-PBS. Standard: Eg6351. Range: 0.312-40 ng/mL