

For Research Use Only

UROM Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87353-3-PBS



Basic Information

Catalog Number: 87353-3-PBS	GenBank Accession Number: NM_001008389.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7369	CloneNo.: 252567C12
Isotype: IgG	UNIPROT ID: P07911-1	
Immunogen Catalog Number: EG5939	Full Name: uromodulin	
	Calculated MW: 70 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of six individual healthy human donors was measured. The Uromodulin concentration of detected samples was determined to be 165.17 ng/mL with a range of 122.41-209.38 ng/mL.



Sandwich ELISA standard curve of MP02973-2, Human Uromodulin Recombinant Matched Antibody Pair - PBS only. 87353-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5939. 87353-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.