

For Research Use Only

Mouse TL1A/TNFSF15 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87367-1-PBS



Basic Information

Catalog Number: 87367-1-PBS	GenBank Accession Number: /	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 326623	CloneNo.: 250983C12
Isotype: IgG	UNIPROT ID: Q5UBV8	
Immunogen Catalog Number: EG5273	Full Name: tumor necrosis factor (ligand) superfamily, member 15	
	Calculated MW: 28 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean TNFSF15 concentration was determined to be 5.70 ng/mL in mouse large intestine extract based on a 4.5 mg/mL extract load.



Serum of six mice was measured. The mouse TNFSF15 concentration of detected samples was determined to be 2.66 ng/mL with a range of 0.87 - 5.86 ng/mL.



Sandwich ELISA standard curve of MP02970-1, Mouse TNFSF15 Recombinant Matched Antibody Pair - PBS only. 87367-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5273. 87367-1-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL