

For Research Use Only

# B3GALT6 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87370-1-PBS



## Basic Information

Catalog Number:

87370-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_080605

GeneID (NCBI):

126792

UNIPROT ID:

Q96L58

Full Name:

UDP-Gal:betaGal beta 1,3-galactosyltransferase polypeptide 6

Calculated MW:

37 kDa

Observed MW:

37-40 kDa

Purification Method:

Protein A purification

CloneNo.:

252785F4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87370-1-RR (B3GALT6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87370-1-PBS in a different storage buffer formulation.