

For Research Use Only

VHL Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87374-2-PBS



Basic Information

Catalog Number: 87374-2-PBS	GenBank Accession Number: NM_000551.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7428	CloneNo.: 251968C9
Isotype: IgG	UNIPROT ID: P40337	
Immunogen Catalog Number: EG6405	Full Name: von Hippel-Lindau tumor suppressor	
	Calculated MW: 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean VHL concentration was determined to be 19.5 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 11.0 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 7.9 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02982-1, VHL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87374-3-PBS. Detection antibody: 87374-2-PBS. Standard: Eg6405. Range: 0.312-40 ng/mL.