

For Research Use Only

VHL Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87374-6-PBS



Basic Information

Catalog Number: 87374-6-PBS	GenBank Accession Number: NM_000551.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7428	CloneNo.: 251968B4
Isotype: IgG	UNIPROT ID: P40337	
Immunogen Catalog Number: EG6405	Full Name: von Hippel-Lindau tumor suppressor	
	Calculated MW: 24 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean VHL concentration was determined to be 37.97 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 17.69 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 4.97 ng/mL in DU145 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02982-2, Human VHL Recombinant Matched Antibody Pair - PBS only. 87374-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6405. 87374-5-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL