

For Research Use Only

CCL4L1 Recombinant monoclonal antibody, PBS Only

Catalog Number:87376-1-PBS



Basic Information

Catalog Number:

87376-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5335

GenBank Accession Number:

N/A

GeneID (NCBI):

9560

UNIPROT ID:

Q8NHW4-1

Full Name:

chemokine (C-C motif) ligand 4-like 1

Calculated MW:

10 kDa

Observed MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

251685H1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PMA, LPS and Brefeldin A treated THP-1 cells were subjected to SDS PAGE followed by western blot with 87376-1-RR (CCL4L1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87376-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87376-1-RR against Human CCL4L1 were performed. The affinity constant is 4.87 nM.