

For Research Use Only

TRIM16 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87377-3-PBS



Basic Information

Catalog Number: 87377-3-PBS	GenBank Accession Number: BC015674	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10626	CloneNo.: 252519D5
Isotype: IgG	UNIPROT ID: O95361	
Immunogen Catalog Number: AG21406	Full Name: tripartite motif-containing 16	
	Calculated MW: 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TRIM16 concentration was determined to be 1654.2 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 47.9 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02988-1, TRIM16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87377-3-PBS. Detection antibody: 87377-2-PBS. Standard: Ag21406. Range: 1.562-200 ng/mL