

For Research Use Only

Human IgM Recombinant monoclonal antibody

Catalog Number: 87380-10-RR



Basic Information

Catalog Number:

87380-10-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6742

GenBank Accession Number:

BC020240

GeneID (NCBI):

3507

Full Name:

immunoglobulin heavy constant mu

Purification Method:

Protein A purification

CloneNo.:

253289A9

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Specificity of the recombinant rabbit anti-human IgM antibody was confirmed by ELISA. The antibody specifically recognizes both human IgM monomer and pentamer, with no cross-reactivity observed against other human immunoglobulins (IgG, IgD, IgE) or IgM from mouse, pig, rat, and cat.