

For Research Use Only

Human IgM Recombinant monoclonal antibody, PBS Only

Catalog Number: 87380-8-PBS



Basic Information

Catalog Number:

87380-8-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6742

GenBank Accession Number:

BC020240

GeneID (NCBI):

3507

Full Name:

immunoglobulin heavy constant mu

Purification Method:

Protein A purification

CloneNo.:

253289F4

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Specificity of the recombinant rabbit anti-human IgM antibody was confirmed by ELISA, demonstrating exclusive reactivity to human IgM monomer with no cross-reactivity to other human immunoglobulins or IgM from mouse, pig, rat, or cat.