

For Research Use Only

ABHD14B Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87389-3-PBS



Basic Information

Catalog Number: 87389-3-PBS	GenBank Accession Number: BC007234	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 84836	CloneNo.: 252778G2
Isotype: IgG	UNIPROT ID: Q96IU4	
Immunogen Catalog Number: AG15097	Full Name: abhydrolase domain containing 14B	
	Calculated MW: 210 aa, 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ABHD14B concentration was determined to be 235.1 ng/mL in L02 cell extract based on a 1.2 mg/mL extract load, 162.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 30.7 ng/mL in U-251 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03004-1, ABHD14B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87389-3-PBS. Detection antibody: 87389-2-PBS. Standard: Ag15097. Range: 0.781-100 ng/mL