

For Research Use Only

ZNHIT1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87390-2-PBS



Basic Information

Catalog Number: 87390-2-PBS	GenBank Accession Number: BC017333	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10467	CloneNo.: 252815B9
Isotype: IgG	UNIPROT ID: O43257	
Immunogen Catalog Number: AG9741	Full Name: zinc finger, HIT type 1	
	Calculated MW: 154 aa, 18 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ZNHIT1 concentration was determined to be 12.1 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 13.9 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 1.3 ng/mL in SKOV-3 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP03003-1, ZNHIT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87390-1-PBS. Detection antibody: 87390-2-PBS. Standard: Ag9741. Range: 0.625-80 ng/mL.