

For Research Use Only

HVEM/TNFRSF14 Recombinant monoclonal antibody

Catalog Number: 87402-1-RR



Basic Information

Catalog Number: 87402-1-RR	GenBank Accession Number: NM_178931.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 230979	CloneNo.: 252718C4
Isotype: IgG	UNIPROT ID: Q80WM9	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: EG3145	Full Name: tumor necrosis factor receptor superfamily, member 14 (herpesvirus entry mediator)	
	Calculated MW: 30 kDa	
	Observed MW: 30-40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse thymus tissue, mouse spleen tissue
Species Specificity: mouse	

Background Information

The mammalian tumor necrosis factor receptor (TNFR) family consists of 10 cell-surface proteins that regulate development and homeostasis of the immune system. Member 14 of TNFR (TNFRSF14, also named as TR2, ATAR, HVEA, HVEM, LIGHTR) was also identified as a mediator of herpesvirus entry into mammalian cells. The cytoplasmic region of TNFRSF14 bound to several members of the TNFR-associated factor (TRAF) family, namely TRAF1, TRAF2, TRAF3, and TRAF5.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87402-1-RR (Tnfrsf14 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87402-1-RR against Mouse Tnfrsf14 were performed. The affinity constant is 30.8 nM.