

For Research Use Only

ST6GAL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87412-1-PBS



Basic Information

Catalog Number: 87412-1-PBS	GenBank Accession Number: BC040009	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6480	CloneNo.: 251718B9
Isotype: IgG	UNIPROT ID: P15907	
Immunogen Catalog Number: EG5012	Full Name: ST6 beta-galactosamide alpha-2,6-sialyltransferase 1	
	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ST6GAL1 concentration was determined to be 4,121.9 pg/mL in Human plasma.



Cytometric bead array standard curve of MP02980-1, ST6GAL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87412-2-PBS. Detection antibody: 87412-1-PBS. Standard: Eg5012. Range: 0.625-80 ng/mL.