

For Research Use Only

ST6GAL1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87412-2-PBS



Basic Information

Catalog Number:

87412-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5012

GenBank Accession Number:

BC040009

GeneID (NCBI):

6480

UNIPROT ID:

P15907

Full Name:

ST6 beta-galactosamide alpha-2,6-sialyltransferase 1

Calculated MW:

47 kDa

Purification Method:

Protein A purification

CloneNo.:

251718E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ST6GAL1 concentration was determined to be 4,121.9 pg/mL in Human plasma.



Cytometric bead array standard curve of MP02980-1, ST6GAL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87412-2-PBS. Detection antibody: 87412-1-PBS. Standard: Eg5012. Range: 0.625-80 ng/mL.