

For Research Use Only

Mouse ITAM Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87419-7-PBS



Basic Information

Catalog Number:	87419-7-PBS	GenBank Accession Number:	Purification Method:
Source:	Rabbit	GeneID (NCBI):	Protein A purification
Isotype:	IgG	16409	CloneNo.:
Immunogen Catalog Number:	EG10463	Full Name:	252750E6
		integrin alpha M	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD11b concentration was determined to be 2,640.3 pg/mL in mouse lung tissue extract based on a 5.0 mg/mL, 24,896.3 pg/mL in mouse spleen tissue extract based on a 6.4 mg/mL extract load and 1,148.8 pg/mL in mouse brain tissue extract based on a 6.3 mg/mL extract load.



Sandwich ELISA standard curve of MP03033-2, Mouse CD11b Recombinant Matched Antibody Pair - PBS only. 87419-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg10463. 87419-7-PBS was HRP conjugated as the detection antibody. Range: 39-2500 pg/mL