

For Research Use Only

HGS Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87420-3-PBS



Basic Information

Catalog Number:

87420-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0589

GenBank Accession Number:

BC003565

GeneID (NCBI):

9146

UNIPROT ID:

O14964

Full Name:

hepatocyte growth factor-regulated
tyrosine kinase substrate

Calculated MW:

86 kDa

Purification Method:

Protein A purification

CloneNo.:

252813A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HGS concentration was determined to be 346.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 327.2 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load, mouse skeletal muscle tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03019-1, HGS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87420-3-PBS. Detection antibody: 87420-2-PBS. Standard: Ag0589. Range: 0.781-100 ng/mL.