

For Research Use Only

# HEATR2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87423-2-PBS



## Basic Information

Catalog Number:

87423-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20080

GenBank Accession Number:

BC047240

GeneID (NCBI):

54919

UNIPROT ID:

Q86Y56

Full Name:

HEAT repeat containing 2

Calculated MW:

855 aa, 94 kDa

Purification Method:

Protein A purification

CloneNo.:

25281488

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean HEATR2 concentration was determined to be 20.3 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load, 7.5 ng/mL in TT cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03032-1, HEATR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87423-3-PBS. Detection antibody: 87423-2-PBS. Standard: Ag20080. Range: 0.781-100 ng/mL