

For Research Use Only

HEATR2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87423-3-PBS



Basic Information

Catalog Number:

87423-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20080

GenBank Accession Number:

BC047240

GeneID (NCBI):

54919

UNIPROT ID:

Q86Y56

Full Name:

HEAT repeat containing 2

Calculated MW:

855 aa, 94 kDa

Purification Method:

Protein A purification

CloneNo.:

252814H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HEATR2 concentration was determined to be 20.3 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load, 7.5 ng/mL in TT cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03032-1, HEATR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87423-3-PBS. Detection antibody: 87423-2-PBS. Standard: Ag20080. Range: 0.781-100 ng/mL