

For Research Use Only

POFUT1/FUT12 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87434-2-PBS



Basic Information

Catalog Number: 87434-2-PBS	GenBank Accession Number: NM_015352.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23509	CloneNo.: 252682C5
Isotype: IgG	UNIPROT ID: Q9H488-1	
Immunogen Catalog Number: EG6687	Full Name: protein O-fucosyltransferase 1	
	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean POFUT1 concentration was determined to be 89.9 ng/mL in CaCO-2 cell extract based on a 1.3 mg/mL extract load, 9.2 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03009-1, POFUT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87434-2-PBS. Detection antibody: 87434-1-PBS. Standard: Eg6687. Range: 0.625-80 ng/mL.