

For Research Use Only

CCL18 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87440-1-PBS



Basic Information

Catalog Number: 87440-1-PBS	GenBank Accession Number: BC069700	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6362	CloneNo.: 251591A4
Isotype: IgG	UNIPROT ID: P55774	
Immunogen Catalog Number: EG5093	Full Name: chemokine (C-C motif) ligand 18 (pulmonary and activation-regulated)	
	Calculated MW: 89 aa, 10 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CCL18 concentration was determined to be 464.4 ng/mL in Human plasma.



Cytometric bead array standard curve of MP03008-1, CCL18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87440-2-PBS. Detection antibody: 87440-1-PBS. Standard: Eg5093. Range: 0.625-40 ng/mL.