

For Research Use Only

FBXO7 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87451-1-PBS



Basic Information

Catalog Number: 87451-1-PBS	GenBank Accession Number: BC008361	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 25793	CloneNo.: 252859B2
Isotype: IgG	UNIPROT ID: Q9Y3I1	
Immunogen Catalog Number: AG1035	Full Name: F-box protein 7	
	Calculated MW: 59 kDa	
	Observed MW: 75 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

FBXO7, also named as FBX7, is a substrate recognition component of a (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. FBXO7 recognizes BIRC2 and DLGAP5. It function in phosphorylation-dependent ubiquitination. The observed molecular weight is consistent with the findings reported in PMID: 26310625 and PMID: 33774278.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87451-1-RR (FBXO7 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87451-1-PBS in a different storage buffer formulation.