

For Research Use Only

NCOA6 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87454-3-PBS



Basic Information

Catalog Number:

87454-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18737

GenBank Accession Number:

BC136272

GeneID (NCBI):

23054

UNIPROT ID:

Q14686

Full Name:

nuclear receptor coactivator 6

Calculated MW:

2063 aa, 219 kDa

Purification Method:

Protein A purification

CloneNo.:

252826E2

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




The mean NCOA6 concentration was determined to be 8.9 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, 12.0 ng/mL in HT-29 cell extract based on a 1.1 mg/mL extract load and 6.4 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP03030-2, Human NCOA6 Recombinant Matched Antibody Pair - PBS only. 87454-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18737. 87454-4-PBS was HRP conjugated as the detection antibody. Range: 19.5-1250 pg/mL.



The mean NCOA6 concentration was determined to be 69.3 ng/mL in CaCO-2 cell extract based on a 1.2 mg/mL extract load, 44.2 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, 13.4 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03030-1, NCOA6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87454-3-PBS. Detection antibody: 87454-2-PBS. Standard: Ag18737. Range: 0.625-80 ng/mL