

For Research Use Only

uPA/Urokinase Recombinant monoclonal antibody, PBS Only

Catalog Number:87458-1-PBS



Basic Information

Catalog Number: 87458-1-PBS	GenBank Accession Number: NM_002658.6	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5328	CloneNo.: 252624A5
Isotype: IgG	UNIPROT ID: P00749-1	
Immunogen Catalog Number: EG7083	Full Name: plasminogen activator, urokinase	
	Calculated MW: 49kDa	
	Observed MW: 49 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

PLAU(Urokinase-type plasminogen activator) is also named as uPA and belongs to the peptidase S1 family.UPA initiates a proteolytic cascade, which degrades extracellular matrix during tissue growth and remodelling.The secreted proform of UPA is activated after binding to a specific high affinity cell surface receptor. Secreted UPA is exposed to rapid inactivation in the pericellular space by plasminogen activator inhibitor-1 (PAI-1), which is stored in the extracellular matrix bound to vitronectin(PMID:9665342).It has 2 isoforms produced by alternative splicing. Defects in PLAU are the cause of Quebec platelet disorder (QPD).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87458-1-RR (PLAU/uPA antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87458-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed PC-3 cells using uPA/Urokinase antibody (87458-1-RR, Clone: 252624A5) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 87458-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87458-1-RR against Human uPA/Urokinase were performed. The affinity constant is 0.549 nM.