

For Research Use Only

Mouse Prss2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87469-2-PBS



Basic Information

Catalog Number: 87469-2-PBS	GenBank Accession Number: NM_009430.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 22072	CloneNo.: 250391H4
Isotype: IgG	UNIPROT ID: P07146	
Immunogen Catalog Number: EG3056	Full Name: protease, serine, 2	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Trypsin 2/PRSS2 concentration was determined to be 8581.8 ng/mL in Mouse pancreas tissue extract based on a 10.6 mg/mL extract load, 1.1 ng/mL in Mouse brain tissue extract based on a 3.2 mg/mL extract load.



Cytometric bead array standard curve of MP03038-1, MOUSE Trypsin 2/PRSS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87469-3-PBS. Detection antibody: 87469-2-PBS. Standard: Eg3056. Range: 0.156-20 ng/mL.