

For Research Use Only

NFASC Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87477-1-PBS



Basic Information

Catalog Number: 87477-1-PBS	GenBank Accession Number: NM_001005388.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23114	CloneNo.: 252657A2
Isotype: IgG	UNIPROT ID: O94856-1	
Immunogen Catalog Number: EG5530	Full Name: neurofascin homolog (chicken)	
	Calculated MW: 150 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Neurofascin concentration was determined to be 32.2 ng/mL in U-87MG cell extract based on a 1.2 mg/mL extract load, THP-1 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03052-1, Neurofascin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87477-2-PBS. Detection antibody: 87477-1-PBS. Standard: Eg5530. Range: 0.625-80 ng/mL.