

For Research Use Only

LY6K Recombinant monoclonal antibody, PBS Only

Catalog Number: 87480-1-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 87480-1-PBS | GenBank Accession Number: NM_017527.3 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 54742 | CloneNo.: 252714A11 |
| Isotype: IgG | UNIPROT ID: Q17RY6-1 | |
| Immunogen Catalog Number: EG5676 | Full Name: lymphocyte antigen 6 complex, locus K | |
| | Calculated MW: 19 kDa | |
| | Observed MW: 18-26 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

The LY6K (Lymphocyte antigen 6K) protein belongs to the LY6/uPAR superfamily and is a class of small glycoproteins that are anchored to the outer side of the cell membrane via glycosylphosphatidylinositol (GPI).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87480-1-RR (LY6K antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87480-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 87480-1-RR (LY6K antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87480-1-PBS in a different storage buffer formulation.