

For Research Use Only

GPR183 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87492-1-PBS



Basic Information

Catalog Number:

87492-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC020752

GeneID (NCBI):

1880

UNIPROT ID:

P32249

Full Name:

G protein-coupled receptor 183

Calculated MW:

361 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

252701B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GPR183 concentration was determined to be 77.7 ng/mL in HL-60 cell extract based on a 1.2 mg/mL extract load, 42.5 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load, 35.1 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03015-1, GPR183 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87492-2-PBS. Detection antibody: 87492-1-PBS. Standard: Ag3035. Range: 1.562-200 ng/mL