

For Research Use Only

PPIL2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87496-1-PBS



Basic Information

Catalog Number: 87496-1-PBS	GenBank Accession Number: NM_148176.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23759	CloneNo.: 252848F8
Isotype: IgG	UNIPROT ID: Q13356-2	
Immunogen Catalog Number: EG6018	Full Name: peptidylprolyl isomerase (cyclophilin)-like 2	
	Calculated MW: 59 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Peptidylprolyl isomerase (cyclophilin)-like 2 (PPIL2), also known as Cyclophilin-like protein Cyp-60, is a U-box-type E3 ubiquitin ligase belonging to the cyclophilin family. PPIL2 is also a target of the JAK2/STAT5 pathway and upregulated in patients with MPNs (PMID: 36092721, 25992900). PPIL2 is a large ~60 kDa (520 amino acids) protein that has a central cyclophilin-like domain and also has an N-terminal U-box domain that exhibits E3 ligase activity. Western detected PPIL2 at an apparent molecular mass of 60-70 kDa (PMID: 40338661, 29352246).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87496-1-RR (PPIL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87496-1-PBS in a different storage buffer formulation.



HEK-293 cells were subjected to SDS PAGE followed by western blot with 87496-1-RR (PPIL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87496-1-PBS in a different storage buffer formulation.