

For Research Use Only

# Mouse PDGF-BB Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87506-2-PBS



## Basic Information

Catalog Number:

87506-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4856

GenBank Accession Number:

AAH23427.1

GeneID (NCBI):

18591

UNIPROT ID:

P31240

Full Name:

platelet derived growth factor, B polypeptide

Calculated MW:

25kd

Purification Method:

Protein A purification

CloneNo.:

252712E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PDGF-BB concentration was determined to be 156.8 ng/mL in mouse serum.



Cytometric bead array standard curve of MP03061-1, MOUSE PDGF-BB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87506-3-PBS. Detection antibody: 87506-2-PBS. Standard: Eg4856. Range: 0.312-40 ng/mL.