

For Research Use Only

IL-17F Recombinant monoclonal antibody, PBS Only

Catalog Number: 87510-1-PBS



Basic Information

Catalog Number: 87510-1-PBS	GenBank Accession Number: NM_052872.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 112744	CloneNo.: 252565B7
Isotype: IgG	UNIPROT ID: Q96PD4	
Immunogen Catalog Number: EG2823	Full Name: interleukin 17F	
	Calculated MW: 18 kDa	
	Observed MW: 18-20 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

The interleukin 17 (IL-17) family of cytokines contains 6 structurally related cytokines, IL-17A, IL-17B, IL-17C, IL-17D, IL-17E and IL-17F. IL-17 family plays crucial roles in host defense against microbial organisms and in the development of inflammatory diseases. IL-17A is a pro-inflammatory cytokine that also has the capacity to promote angiogenesis and osteoclastogenesis. IL-17F shares the highest homology with IL-17A and signals via a receptor composed by the IL-17RA and IL-17RC subunits. IL-17A and IL-17F can form IL-17A/A or IL-17F/F homodimers, IL-17A/F heterodimers are also formed. IL-17A and IL-17F, produced by the Th17 CD4(+) T cell lineage, have been linked to a variety of inflammatory and autoimmune conditions. IL-17F levels are elevated in sera and lesional psoriatic skin compared to non-lesional tissue. IL-17F also has been implicated in the development of neutrophilic airway inflammation.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87510-1-RR (IL17F antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87510-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87510-1-RR against Human IL-17F were performed. The affinity constant is 43.1 nM.