

For Research Use Only

Mouse CD318/CDCP1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87514-2-PBS



Basic Information

Catalog Number: 87514-2-PBS	GenBank Accession Number: NM_133974.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 109332	CloneNo.: 251667G9
Isotype: IgG	UNIPROT ID: Q5U462	
Immunogen Catalog Number: EG6199	Full Name: CUB domain containing protein 1	
	Calculated MW: 93 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CDCP1/CD318 concentration was determined to be 13.9 ng/mL in mouse skin tissue extract based on a 3.9 mg/mL extract load, mouse skeletal muscle tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03057-1, MOUSE CDCP1/CD318 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87514-3-PBS. Detection antibody: 87514-2-PBS. Standard: Eg6199. Range: 0.312-40 ng/mL.